


# Ambulance Pipeline Systems

## Assembled - Single



| Gas    | Part No. | Price Each |
|--------|----------|------------|
| O2     | 116191A  | £ 75.00    |
| O2/N2O | 116192A  | £ 75.00    |

Other gases available on request

## Single Ambulance Unassembled Pipeline

Comprises:

- 1 x Ambulance Surface Terminal Unit
- Up to 7m Anti-static Hose
- 1 x Remote Probe
- 6 x Pipe Clips & Screws
- 2 x Hose Ferrules
- 2 x Ferrule Sleeves

| Gas    | Part No. | Price Each |
|--------|----------|------------|
| O2     | 116191   | £ 67.50    |
| O2/N2O | 116192   | £ 67.50    |

Other gases available on request

## Assembled - Twin



| Gas    | Part No. | Price Each |
|--------|----------|------------|
| O2     | 116201A  | £ 130.75   |
| O2/N2O | 116202A  | £ 130.75   |

Other gases available on request

## Twin Ambulance Unassembled Pipeline

Comprises:

- 2 x Ambulance Surface Terminal Units
- Up to 7m Anti-static Hose
- 1 x Remote Probe
- 1 x 'T' Piece
- 6 x Pipe Clips & Screws
- 6 x Hose Ferrules
- 6 x Ferrule Sleeves

| Gas    | Part No. | Price Each |
|--------|----------|------------|
| O2     | 116201   | £ 111.25   |
| O2/N2O | 116202   | £ 111.25   |

Other gases available on request