



# Cylinder Accessories

## Cylinder Trolley



Description	Part No.	Price Each
'E' Size (24 cu ft)	TH8803	£ 99.75
'F' Size (48 cu ft)	TH8802	£ 120.00
'G' Size (120 cu ft)	TH8801	£ 141.00
'J' Size (240 cu ft)	TH8800	£ 215.00


## Cylinder Rings & Bracket




Description	Part No.	Price Each
A. Cylinder Wall Bracket only	TH8822	£ 29.50
B. 'C' Size	TH8823	£ 29.00
B. 'D/E' Size	TH8824	£ 29.00
B. 'F' Size	TH8825	£ 29.00
B. 'G' Size	TH8826	£ 29.00
B. 'J' Size	TH8827	£ 29.00

## Medical Cylinder Chart

	C	D	E	F	G	J
Oxygen	Pin Index	Pin Index	Pin Index	Bull Nose	Bull Nose	Pin Index
Nitrous Oxide	Pin Index	Pin Index	Pin Index	Handwheel	Handwheel	Handwheel
Entonox		Pin Index		Pin Index	Pin Index	
Air			Pin Index	Bull Nose	Bull Nose	Pin Index



Pin Index Valve



Side Spindle Pin Index Valve



Handwheel Valve



Bullnose Valve