





Electric Suction

Standard suction units are fitted with autoclavable polysulphone jar(s) - 1800ml capacity + safety overflow valve.


Disposable bacterial filter is fitted as standard - Hydrophobic filter can be fitted on request.

Portable Electric Suction - 240V AC



Description	Type	Part No.	Price Each
A. 240V Portable High Suction. Handle + Single Jar	Classic	FSE 500-110019	£ 490.00
	Premier	FSE 500-111019	£ 535.00
B. 240V Portable High Suction. Handle + Two Jars	Classic	FSE 550-110021	£ 550.00
	Premier	FSE 550-111021	£ 605.00
C. 240V Mobile High Suction Handle + Two Jars	Classic	FSE 570-120021	£ 645.00
	Premier	FSE 570-121021	£ 695.00
D. 240V Portable Low Suction. Handle + Single Jar	Premier	FSE 600-111033	£ 595.00
E. 240V Portable Intermittant. Low Suction. Handle + Single Jar	Premier	FSE 650-111033	£ 890.00

Portable Electric Suction - 12V DC/240V AC



Description	Part No.	Price Each
Portable Emergency Battery/ Mains High Suction. Single Jar	FSE 450-110013	£ 710.00