

AGSS Flow Rate Test Set

SPECIFICATION SHEET

For Checking AGSS Terminal Units to British Standards BS 6384:1987.



Flow Gauge Features

- Large flow capacity for use
- High visibility
- Instantaneous in operation.

Kit Features

- Flow gauge with AGSS hose for connectivity to outlet point for allowing basic flow rate test to BS 6834:1987
- 1, 2 & 4 kPa Orifices supplied
- In compliance with HTM02-01 and relevant European Standards and all preceding standards
- Supplied in a purpose built lightweight carry case
- MEC Medical Ltd provide an annual calibration service to ensure equipment is maintained to the highest standard.

Specification

Weight	2.2Kg
Length	38cm
Width	32cm
Depth	10cm
Recalibration Ref:	RECALDIS



Supplied by