

Emergency Kit For Pressure Gases

SPECIFICATION SHEET

The Emergency kit is used with an Area Valve Service Unit (AVSU) as a back-up in case of an emergency, should the conventional supply system fail, or for carrying out routine maintenance work.

One kit will serve two cylinders either free standing or in trolleys, enabling a quick change over with the minimum of delay.

All components are gas specific to BS 5682:1998.

Features

- One kit will serve two cylinders
- 2 x Two Stage Bullnose or Pin Index Regulator with right angle adaptor
- 2 x Splitter Hose Assembly
- 1 x Splitter Unit
- 1 x Pendant Hose Assembly
- 1 Emergency Kit Box.

ONLY INCLUDED IN UNIVERSAL KIT:

- 2 x Bull Nose—Pin Index Adaptors
- Combination Spanner
- Cylinder Key.



- Gas specific available for O₂, Medical Air, Surgical Air, Nitrous Oxide and Entonox.
- Bull Nose, Pin Index & Nut and Stem Cylinder Connexions.
- Regulators Pre-set to 4Bar, 7Bar for surgical Air.
- Universal Kit available for cylinder connection to both Bull Nose and Pin Index Cylinders, with the extra inclusion of Bull Nose—Pin Index Adaptors, Cylinder Spanner and Pin Index Key.

Specification

Weight	11.90Kg
Height	45cm
Width	55cm
Depth	20cm



Supplied by