

## F / HX Cylinder Trolley

SPECIFICATION SHEET

**MEC cylinder trolleys are designed to be ergonomically constructed for ease of use.**

The red securing straps are designed to remove all cylinder noise and be a clear and visible indicator of the cylinder being secure.

All of our new trolleys are painted in a non-slip paint which has anti-bacterial properties and is extremely robust. The large wheels aid in the manoeuvring of the trolley with its counter balance properties no effort is required to move a cylinder and its stability in the vertical position is exceptional.


### Specification

Weight: 8.5Kg

Height: 106cms

Width: 36cms

Depth: 32cms



Supplied by