

Flowmeter Range

SPECIFICATION SHEET

MEC Flowmeters are designed and built in the UK to the highest standards with high and low flow options.

MEC Medical Ltd Flowmeter range fully complies with the requirements of ISO15002:2000 Flow metering devices for the connection to terminal units of medical gas pipeline systems. MEC's Flowmeters are calibrated to operate from an inlet pressure of 4 Bar (approx. 60psi) CE marked all MEC Flowmeters carry a 1 year function warranty.

Features

- Anodised Aluminium Body
- Acetal Spindle
- Polycarbonate Inner and Outer Tube
- Nitrile O Rings
- White Acetyl Plastic Tubing Nipple
- Single and Double variants available
- 3/8" BSP Nut & Liner, Direct Probe or Rail mounting options
- High and Low Flows for Oxygen and Medical Air.



Supplied by