


High Performance Emergency Foot/Hand Operated Suction Unit

SPECIFICATION SHEET

MEC Emergency Foot/Hand Operated Suction Units are designed and built to the highest standards for use without the need for any electrical power supply.



FSF 100-110000

Features

- High performance emergency foot operated suction unit
- Sturdy, well-balanced and fully portable with carrying handle
- Superior suction performance without the need for any electrical power supply
- Easy and reliable operation
- Includes a 1800ml Autoclavable receiver jar with an overflow shut-off valve
- The unit includes a high vacuum gauge and vacuum control knob
- Double action piston pump offers enhanced suction performance.

Specification

Maximum Vacuum	-94kPa (708 mmHg)
Collection Jars	1800ml Autoclavable receiver jar with overflow protection, x1 (graduated to 1500ml)
Dimensions	34cm (L) x 15cm (H) x 16cm (W)
Weight	Approx. 4kg



Supplied by