

MRI Compatible Dial Flowmeter / Regulators

SPECIFICATION SHEET

MEC are pleased to announce a new range of Dial Flow Regulators which are non-magnetic and MRI compatible to add to our portfolio. MEC Dial Flowmeter/Regulators are designed and built to the highest standards with clear display of the flow rate and multiple flow therapy output.

Features

- Non-magnetic for use in MRI
- Oxygen regulator with Pin Index, Bullnose or CGA input
- Multiple flow therapy output 0, ½, 1, 2, 3, 4, 6, 8, 10, 15 and MAX 25 litres per minute
- Inlet pressure of 200 bar maximum
- Barbed fir tree style connector
- Designed so that the contents gauge and flow rate display window can both be easily seen, at the same time the therapy outlet and main 4 bar oxygen outlet can be accessed so that equipment may be attached.


Bullnose Regulator
+ Flowmeter



Pin Index Regulator +
Flowmeter



USA Type Regulator +
Flowmeter



Supplied by