

Male NIST Anti Confusion Probes

SPECIFICATION SHEET

Male NIST Anti Confusion Probe features

- Set of six or nine Stainless Steel gas specific NIST probes
- Gases include O2, O2/N2O, N2O, Air-4 Bar, Air-7 Bar, Vacuum
- Nine includes CO2, N, HEO2 (Carbon Dioxide, Nitrogen, Heliox)
- Calibrated and manufactured to current industry standards
- In compliance with HTM02-01 and relevant European Standards and all preceding standards.



Male NIST Anti Confusion Probes

Description	Part Number
Male NIST Anti Confusion Probes x6	20170
Male NIST Anti Confusion Probes x9	20170-9

Specification

Weight:	1.75Kg
Height:	23cm
Depth:	7cm
Width:	27cm
Recalibration Ref:	RECALACN



Supplied by