



# Male and Female Ceiling NISTs

SPECIFICATION SHEET

**Features**

- 150mm Copper Stub, 15mm diameter (Other lengths available on request)
- Gas Specific engraving
- Female NIST includes non-return valve (Excludes Vacuum and AGSS)
- Medical degreased copper tube
- MRI Compatible
- Fully conforms to EN ISO 15908:2010



**Ceiling NIST - Female**

Gas	Part No.
O2	3601
O2/N2O	3602
N2O	3603
Air-7	3604
Air-4	3605
Vacuum	3606
Co2	3607
Nitrogen	3608
AGSS	231724

**Ceiling NIST - Male**

Gas	Part No.
O2	3601N
O2/N2O	3602N
N2O	3603N
Air-7	3604N
Air-4	3605N
Vacuum	3606N
Co2	3607N
Nitrogen	3608N
AGSS	231724N



Supplied by