

Mobile Vacuum Theatre High Suction Unit

SPECIFICATION SHEET

MEC Mobile Electric Suction Controllers are designed and built to the highest standards. Running on 240V AC, supplied with two jars and fitted to castors.

Standard suction units are fitted with autoclavable polysulphone jar(s) – 1800ml capacity and safety overflow valve. Disposable bacterial filter is fitted as standard – Hydrophobic filter can be fitted on request.

Features

- Mobile electric high vacuum suction unit with two lockable castors
- Suitable for use in the operating theatre.

Classic Unit

Also available as
Premier Unit



Specification

Weight	17.5Kg
Height	88cms
Depth	45cms
Width	32cms
Suction Capacity	45 l/min
Max Vacuum	-84kPa (630mmHg)
Filter	Disposable bacterial filter (BIF 001)



Supplied by