

Mobile Electric High Suction – On Trolley

SPECIFICATION SHEET

MEC Mobile Electric Suction Controllers are designed and built to the highest standards. Running on 240V AC, supplied with two jars and fitted to a mobile trolley.

Standard suction units are fitted with autoclavable polysulphone jars – 1800ml capacity and safety overflow valve. Disposable bacterial filter is fitted as standard – Hydrophobic filter can be fitted on request.

Features

- Mobile electric high vacuum suction unit, mounted on a trolley
- Suitable for use in minor procedures.

Classic Model



Also suitable for Premier Model



Specification

Weight	7.2Kg
Height	295mm
Depth	230mm
Width	560mm
Suction Capacity	26 l/min
Power Consumption	75 Watts
Max Vacuum	-83kPa (623mmHg)
Filter	Disposable bacterial filter (BIF 001)



Supplied by