

Pendant Hoses

SPECIFICATION SHEET

MEC supplies a complete range of Medical Gas Hoses which are fully compliant to current European Standards for manufacturing.



FPH1 O2 Pendant Hose Assembly shown
(includes cover, not shown)

MEC Medical Gas Hoses are part numbered, Lot numbered and date stamped for full traceability and include a replacement date. CE marked and insured, the highest build quality is assured with all fittings manufactured in-house. We have the complete range in all lengths with custom lengths also available.

- Gas specific coloured Anti-Static Hose
- MRI Compatible
- Anti-Microbial
- Phthalate Free
- Gas Specific NIST connectors
- Full range of hose lengths available manufactured to the customer's requirement

1m Pendant Hose Assembly:

Gas	Part No.
O2	FPH1
O2/N2O	FPH2
N2O	FPH3
Air-7	FPH4
Air-4	FPH5
Vacuum	FPH6
AGSS	FPHAGSS

Alternative hose length may be specified at time of order.



Supplied by: