

Pin Index Dial Regulator with Flowmeter

SPECIFICATION SHEET

**MEC Medical Pin Index Regulator combined with Flowmeter
available in O₂ or Air-4 gases. Many other configurations available.**

The MEC range of Pin Index Medical Regulators and Flowmeters are CE approved and are manufactured under strict BS EN ISO 13485 quality design and management systems. Designed and manufactured to fully comply with all current standards.

Features

- Nickel plated Pin Index inlet connection
- White or Black acetal plastic tubing nipple on flowmeter
- Acetal spindle on flowmeter
- Robust brass body construction for prolonged service life
- Serial number and replacement date, clearly marked on rear of regulator.
- Multiple flow therapy output 0, ½, 1, 2, 3, 4, 6, 8, 10, 15 and MAX 15 litres per minute

4918511 shown



Specification

Inlet Pressure	200 bar +/-50 bar maximum
Delivery Pressure	0, ½, 1, 2, 3, 4, 6, 8, 10, 15 and MAX 15 litres per minute
Maximum Flow Rate	-15 l/min
Available for a range of medical gases	Oxygen Air-4



Supplied by