

Portable Electric High Suction – Two Jars

SPECIFICATION SHEET

**MEC Electric Suction Controllers are designed and built to the highest standards.
Running on 240V AC and supplied with two jars and a handle.**

Standard suction units are fitted with autoclavable polysulphone jars – 1800ml capacity and safety overflow valve. Disposable bacterial filter is fitted as standard – Hydrophobic filter can be fitted on request.

Features

- Portable electric high vacuum suction unit
- Two collection jars
- Fitted with a change-over valve that enables the operator to switch from one jar to the other without interrupting suctioning.



Premier Model



Classic Model

Specification

Weight	7.2Kg
Height	295mm
Depth	230mm
Width	560mm
Suction Capacity	26 l/min
Power Consumption	75 Watts
Max Vacuum	-83kPa (623mmHg)
Filter	Disposable bacterial filter (BIF 001)



Supplied by