

## Portable Electric High Suction – Single Jar

### SPECIFICATION SHEET

**MEC Electric Suction Controllers are designed and built to the highest standards.  
Running on 12V/240V AC, supplied with a single jar and handle.**

Standard suction unit is fitted with autoclavable polysulphone jar – 1800ml capacity and safety overflow valve. Disposable bacterial filter is fitted as standard – Hydrophobic filter can be fitted on request.

#### Features

- Portable emergency Battery/Mains high vacuum suction unit
- One collection jar
- Can be operated from the built-in battery, a 220V AC or 12V DC source.



#### Specification

Weight	7.1Kg
Height	295mm
Depth	230mm
Width	435mm
Suction Capacity	26 l/min
Power Consumption	75 Watts
Max Vacuum	-80kPa (600mmHg)
Filter	Disposable bacterial filter (BIF 001)



Supplied by