

## Portable Electric Low Suction – Single Jar

### SPECIFICATION SHEET

**MEC Electric Suction Controllers are designed and built to the highest standards.  
Running on 240V AC and supplied with a single jar and a handle.**

Standard suction unit is fitted with autoclavable polysulphone jar – 1800ml capacity and safety overflow valve.

Disposable bacterial filter is fitted as standard – Hydrophobic filter can be fitted on request.

#### Features

- Portable electric low vacuum suction unit
- One collection jar
- Suitable for abdominal and pleural drainage.



#### Specification

Weight	5Kg
Height	295mm
Depth	230mm
Width	435mm
Suction Capacity	4 l/min
Power Consumption	4 Watts
Max Vacuum	-20kPa (150mmHg)
Filter	Disposable bacterial filter (BIF 001)

Supplied only as Premier Model



Supplied by