

Portable Intermittent Low Suction

SPECIFICATION SHEET

**MEC Electric Suction Controllers are designed and built to the highest standards.
Running on 240V AC and supplied with a single jar and a handle.**

Standard suction unit is fitted with autoclavable polysulphone jar – 1800ml capacity and safety overflow valve.
Disposable bacterial filter is fitted as standard – Hydrophobic filter can be fitted on request.

Features

- Portable electric low vacuum suction unit
- One collection jar
- Features continuous and intermittent suctioning and is suitable for abdominal and pleural drainage.



Specification

Weight	6Kg
Height	295mm
Depth	230mm
Width	435mm
Suction Capacity	4 l/min
Power Consumption	4 Watts
Max Vacuum	-20kPa (150mmHg)
Filter	Disposable bacterial filter (BIF 001)

Supplied only as Premier Model



Supplied by