

# MEC

MEDICAL LIMITED

## Purair 500

### SPECIFICATION SHEET

MEC Medical Ltd provides a range of Anaesthetic Gas Scavenging Systems for a variety of applications. Each system differs in its extraction rate with flow capacity matched to the needs of a specific environment.



- Portable AGSS system primarily designed for midwifery applications, can be used in other scenario where gas scavenging is required
- Eliminates the risk of exposure to waste anaesthetic gases, especially when working in confined spaces
- Extraction capacity is 500 litres per minute
- Simple to operate with robust design and lightweight carry case
- Carry Case includes wheels and telescopic handle.

#### Specification

Weight	13Kg
Height	59cm
Depth	26cm
Width	35cm

CE

Supplied by:

# MEC

MEDICAL LIMITED

MEC Medical Limited, Unit 3, Trust Industrial Estate,  
Wilbury Way, Hitchin, Hertfordshire. SG4 0UZ. United Kingdom

T: +44 (0) 1462 436396 F: +44 (0) 1462 432411  
sales@mecmedical.com www.mecmedical.com