

# Single Stage Regulator + Extension Bar

## SPECIFICATION SHEET

**MEC Medical Single Stage Regulators in a wide range of configurations and gases with International inlet connections available.**

The MEC range of Single Stage Medical Regulators are CE approved and are manufactured under strict BS EN ISO 13485 quality design and management systems at our UK production facility. Designed and manufactured to fully comply with BS EN ISO 10524-1 and The Medical Devices Directive 93/42/EEC.

### Features

- Bullnose, Pin Index, Nut & Stem and International inlet connections available
- Nickel plated inlet and outlet connections at either bottom or side entry
- Robust brass body construction for prolonged service life
- Serial number and replacement date, clearly marked on rear of regulator.



6080FEXT Shown

### Specification

Inlet Pressure	200 bar +/-50 bar maximum
Delivery Pressure	4 bar pre-set
Maximum Flow Rate	-50 l/min
Available for a range of medical gases	Oxygen Air Nitrous Oxide Entonox



Supplied by