

Standard Particulate & Purging Kit

SPECIFICATION SHEET



- Designed for testing at outlet/terminal unit
- Polycarbonate filter holder 47mm
- 75/150Lpm Flow rate jets
- Hydrophobic Membrane 47mm filter papers, 0.45(µm) Pore size
- In compliance with HTM02-01 and relevant European Standards and all preceding standards.

Supplied in a purpose built lightweight carry case.

MEC Medical Ltd provide an annual calibration service to ensure equipment is maintained to the highest standard.

Specification

Weight	2.00Kg
Length	38cm
Width	32cm
Depth	10cm
Recalibration Ref:	RECALPURG



Supplied by