

Anti Confusion Probes

SPECIFICATION SHEET



Anti Confusion Probe features

- Set of six stainless steel gas specific probes
- Gases include O₂, O₂/N₂O, N₂O, Air-4 Bar, Air-7 Bar, Vacuum
- Calibrated to BS5682:19998
- In compliance with HTM02-01 and relevant European Standards and all preceding standards.

Anti Confusion Probes

Description	Part Number
Anti Confusion Probes	20160

Specification

Weight:	1.25Kg
Height:	23cm
Depth:	7cm
Width:	27cm
Recalibration Ref:	RECALCON



Supplied by