

## Terminal Units

### SPECIFICATION SHEET

**MEC Terminal Units are designed and built in the UK to the highest standards with integral non-interchangeability.**



MEC Medical Ltd Terminal Units are designed to fully comply with the requirements of European Standards BS EN ISO 9170-1 – 2008. MEC Terminal Units can be either surface flush or trunking mounted. Features integral to the units include non-interchangeability between probes, socket assemblies and terminal blocks.

Quick release of the probe with single handed operation. Incorporates a secondary valve, which facilitates ease of maintenance without interruption of the gas supply to the unit.

#### Features

- Gas specific sockets with non-interchangeability
- Complete with fixing screws to secure both First Fix and Box to base plate
- First and Second Fix Supplied
- Fully compliant with BS EN ISO 9170-1 – 2008
- All gases available including AGSS.

#### Additional Features

- Mounting styles including Flush, Surface, Trunking and Boom types
- First fix can be supplied with copper stub, Female or Male NIST fittings and Ambulance Hose Tail.



Supplied by