


Twin Adaptor On Straight Bar with British Standard Probe

SPECIFICATION SHEET

MEC Twin Adaptors On Straight Bar with British Standard Probes are made to the highest standards with a choice of all pressure gases available.



201901B Shown

MEC Medical Ltd Twin Adaptors on Straight Bar with British Standard Probes comply to all current safety standards and can be supplied in all gas configurations.

Features

- Non return gas specific Schrader valves allow the split of gas from a single outlet point
- Available with a Direct British Standard BS Probe
- Parallel Configuration
- Lightweight design to connect to a British Standard Outlet point.

Twin Adaptor On Straight Bar with British Standard Probe:

Gas	Part No.
O2	201901B
O2/N2O	201902B
N2O	201903B
Air-7	201904B
Air-4	201905B
Vacuum	201906B

Specification

Weight	0.5Kg
Width	10cm
Depth	13cm
Height	4cm



Supplied by: