

## Twin Adaptor with 3/8" BSP Nut and Liner

### SPECIFICATION SHEET

MEC Twin Adaptors with 3/8" BSP Nut and Liners are made to the highest standards with a choice of all pressure gases available.



201901 Shown

MEC Medical Ltd Twin Adaptors with 3/8" BSP Nut and Liners comply to all current safety standards and can be supplied in all gas configurations.

#### Features

- Non return gas specific Schrader valves allow the split of gas from a single outlet point
- Available with a 3/8" Nut and Liner
- 'Y' Configuration
- Lightweight design to connect to a British Standard Outlet point.

#### Twin Adaptor with 3/8" BSP Nut and Liner:

| Gas    | Part No. |
|--------|----------|
| O2     | 201901   |
| O2/N2O | 201902   |
| N2O    | 201903   |
| Air-7  | 201904   |
| Air-4  | 201905   |
| Vacuum | 201906   |

#### Specification

|        |       |
|--------|-------|
| Weight | 0.5Kg |
| Width  | 13cm  |
| Depth  | 13cm  |
| Height | 4cm   |



Supplied by: