

Twin Adaptor with British Standard Probe

SPECIFICATION SHEET

MEC Twin Adaptors with British Standard Probes are made to the highest standards with a choice of all pressure gases available.



201901P Shown

MEC Medical Ltd Twin adaptors with British Standard probes comply to all current safety standards and can be supplied in all gas configurations.

Features

- Non return gas specific Schrader valves allow the split of gas from a single outlet point
- Available with a Direct British Standard BS Probe
- 'Y' Configuration
- Lightweight design to connect to a British Standard Outlet point.

Twin Adaptor with British Standard Probe:

Gas	Part No.
O2	201901P
O2/N2O	201902P
N2O	201903P
Air-7	201904P
Air-4	201905P
Vacuum	201906P

Specification

Weight	0.5Kg
Width	13cm
Depth	13cm
Height	4cm



Supplied by: